



Mini Portfolio Folder

ex2-1253 / A4 / PATTERN 1 / 1 of 2



50 mm reference

Print the matching paper size at 100% / actual size.
Measure this line before cutting, folding or tracing.

Mini Portfolio Folder

ex2-1253 / A4 / INSTRUCTIONS / 2 of 2

A small two-panel folder with an unglued spine and two folded lower pockets for paper samples.

DIMENSIONS

Two main panels: 75 x 90 mm; spine: 5 mm.

Lower pockets: 30 mm high.

Outside joining tabs: 10 mm.

MAKE A PLAIN-PAPER TRIAL

- 01 Cut the complete outline and the small separation at the bottom of the spine.
- 02 Score both spine lines, the lower pocket folds and the two outer tab folds.
- 03 Fold both lower pockets upward. Fold their outer tabs around the back of the corresponding main panel.
- 04 Glue only those outer tabs. Keep the pocket mouths and spine clear.
- 05 Insert two small paper samples, close the folder and check that the pockets do not catch at the spine.


MATERIALS

Ordinary paper for a first trial; scissors; ruler; a blunt scoring tool; a small amount of paper adhesive.

BEFORE YOU USE IT

Miniature paper-sample folder, not an archival portfolio or commercial sheet-size fit.

Solid boundary = cut. Long dashes = fold/score. Dots = positioning guide. Pale grey = joining area. Keep tabs attached. Adult help is needed for internal cuts.



Print the matching paper size at 100% / actual size.
Measure this line before cutting, folding or tracing.