


Seed-sharing packet with collection fields

Nominal panel 60 x 90 mm / Solid = cut, dashed = fold / Trial-fit your contents



See the second page for folding, gluing and trial-fit instructions.

Seed-sharing packet with collection fields

Construction notes / Use with the matching template page

01 / Prepare

Print at actual size and measure the 50 mm line. Make the first packet with a spare sheet. Complete Seed, Variety, Collected and Shared by while the sheet is flat.

02 / Cut

Cut the solid outer silhouette. Keep both narrow side tabs, the lower front panel and the top closure flap attached. Do not cut the dashed fold lines.

03 / Dry-fold

Crease the four dashed lines. Fold the two 10 mm side tabs inward, then bring the lower 80 mm front panel up against them. Check the printed face and opening before applying adhesive.

04 / Join

Apply a thin line of paper adhesive where each side tab meets the front panel. Press the two side joins and let them dry. Leave the top opening and 28 mm closure flap free.

05 / Check the pocket

The back panel is 60 x 90 mm; folds reduce the usable space. Check the bottom corners and side joins for gaps. Trial the packet over a tray with a few seeds of the intended size.

06 / Fill and close

Use dry seeds and leave room to fold the top flap. Close it over the opening and secure it. Add a small inner paper fold if seeds escape during the trial; label it with the same seed details.



50 mm scale check | Measure after printing